# AL-SENSE™ kit

## Amniotic Leak Self-Test Detector

**Carefully read the entire insert before using this product.**

### The Amniotic Leak Self-Test Detector Kit

The AL-SENSE™ detector is an in vitro self-test, easy-to-use home test kit that allows you to detect probable leakage of amniotic fluid whenever your physician suspects you may be in a higher risk than normal of amniotic fluid leak. Risk of amniotic fluid leak is higher than normal during 48 hours after amniocentesis tests, or when encountering **unidentified wetness**, or among women with **history of premature membrane rupture**. The AL-SENSE™ changes colors when it comes in contact with fluid of pH levels greater than or equal to 6.5 units. Amniotic fluid pH levels are greater than 6.5.

### Principle of Operation

Simply attach the special AL-SENSE™ to your underwear and get on with your day. Wait until you feel a fluid leak, remove the panty liner from the underwear, and check for a stain (any intensity, shape, size or location) on the yellow background. If the panty liner stains light blue or green, amniotic fluid may be leaking from your uterus, or you may have a vaginal infection.

You should consult your physician or go to the hospital without delay. If the panty liner remains yellow, the fluid leak is probably urine, which is common during pregnancy. Any minimal amounts of amniotic fluid leakage which can be sensed as wetness by the user will create a visible blue or green stain. For further diagnosis and medical care, report the results to your physician.

The AL-SENSE™ is recommended for screening either high-risk or normal pregnancies. Leaking amniotic fluid increases the risk of infection to both mother and fetus and early detection can help to:

- Prevent complications or premature birth.
- Identify a possible membrane rupture, especially after amniocentesis.
- Confirm that your water has broken so that you arrive at the hospital in time for a safe delivery.

### Performance Data:

A total of 412 pregnant women were enrolled in 2 different clinical studies of the AL-SENSE Test Kit (103 and 309 pregnant women). The AL-SENSE, as a continuance test, demonstrated a sensitivity ranging from 95.7% to 100% (1, 2). The clinical diagnosis of amniotic leak that was done in the hospitals reached a sensitivity of 97.6% (2).

1) Bornstein J, Geva A, Salt I, Fait V, Schoenfeld A, Shoham HK, Sobel J. Noninvasive diagnosis of premature ruptured amniotic membranes using a novel polymer. Am J Perinatol. 2006 Aug; 23(6):351-4.

2) Bornstein J, Ohel G, Sorokin Y, Reape KZ, Shnaider O, Kessary-Shoham H, Ophir E. Effectiveness of a novel home-based testing device for the detection of rupture of membranes. Am J Perinatol. 2009 Jan;26(1):45-50.

**The AL-SENSE™ gives you a little extra peace of mind — when you need it most.**

### Indications for Use

This product is intended to detect leaking amniotic fluid and identify cause of wetness during pregnancy. For further diagnosis and medical care, report or show test results to your physician.

**If you suspect you have a vaginal infection, consult your Physician.**

### Limitation of the device

Antibiotic therapy or vaginal infections (such as bacterial vaginosis or trichomoniasis) can lead to an elevated vaginal pH level, which may result in a false-positive test result for the presence of amniotic fluid. The possibility of a vaginal infection may be evaluated using standard diagnostic procedures.

### Package Contents

The AL-SENSE™ contains 1 test panty liners, and instructions for use.

Keep the AL-SENSE™ sealed and use it immediately after opening.

### Conditions for Use

To ensure valid results, do not use the AL-SENSE™ monitoring device if the answer to any of the following questions is YES:

- (1) is it less than 12 hours after sexual intercourse?
- (2) is it less than 12 hours after vaginal douching?

During bleeding results reading may need professional assistance. Report all bleeding situations to your physician.

Disregard positive result and retest if there is a possibility of wetness/wetness remains from shower, bath etc, as these may interfere with the results.

In case of sensitization or irritation of skin, at the application site, stop use and contact your physician.

### Directions for Use

1. Tap water can interfere with the test — wash and dry your hands.
2. Unwrap the AL-SENSE™ panty liner. The panty liner contains a pale yellow indicator.
3. Attach the panty liner to your underwear (Figure 1). The pale yellow indicator strip should sit directly against your vagina. Wear the panty liner as you would a normal panty liner till you feel wetness (do not exceed 12 hours).


Figure 1

4. At the first sense of vaginal wetness, or when replacing with a new panty liner, remove the panty liner from your underwear, and check for a stain on the panty liner's yellow background.
5. If you feel wetness and the panty liner is completely yellow, after using it, consider using another test.

### Reading your Results

Immediately after removing the panty liner look for a stain on the panty liner (Figure 2). Result reading should be done in a well illuminated place.

Probable amniotic fluid leak

Urine leak

Figure 2

If the panty liner is stained blue or green (any intensity), in any size, shape or location, on the panty liner, the fluid leak is probably amniotic fluid. The stain can look cloudy and can be more intense along part of the yellow strip edge. In any such case, consult your physician or go to the hospital without delay.

**Important:** The indicator strip may also turn blue or green if an infectious vaginal disease is detected. If the indicator strip turns blue or green, always consult your physician.

If the panty liner did not change color (remained yellow), the vaginal wetness was probably not caused by amniotic fluid. If you continue to experience vaginal wetness, apply a new panty liner and repeat the test or contact your healthcare provider.

### Storage

Store in a dry place at room temperature. Keep the AL-SENSE™ sealed until use to avoid contamination. Please note the expiration date stated on the AL-SENSE™. If in doubt, consult with the dispensing pharmacist.

Keep out of reach of children.

## User Assistance Information and Distribution Details

### Info of distributor

### EU Authorized Representative

Medical Device & QA Services  
76 Stockport Road., Timperley Cheshire, WA15 7 SN, United Kingdom  
Fax: +44 161 903 9787, e-mail: bjt@tremaine.fsnet.co.uk

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CS Common Sense Ltd.